

| Modell                                            | USL        |            |            |            |            |            |
|---------------------------------------------------|------------|------------|------------|------------|------------|------------|
| Speedy 4all Ergo AF5, kurz                        |            |            | 37-38      | 39-43      | 44-47      | 48-50      |
| Speedy 4all Ergo AF5, lang                        |            |            | 37-38      | 39-43      | 44-48      | 49-50      |
| Speedy 4all AF4-4, lang und kurz, Kniewinkel 95°  |            |            | 37-41      | 42-46      | 47-50      |            |
| Speedy 4all AF4-4, lang und kurz, Kniewinkel 105° |            |            | 37-45      | 46-50      |            |            |
| Speedy A1 und A2 kurz, Kniewinkel 95°             |            |            | 33-38      | 39-43      | 44-47      | 48-50      |
| Speedy A1 und A2 kurz, Kniewinkel 105°            |            |            | 33-39      | 40-44      | 45-48      | 49-50      |
| Speedy A1 und A2 lang, Kniewinkel 95°             |            |            | 35-39      | 40-44      | 45-49      | 50-52      |
| Speedy A1 und A2 lang, Kniewinkel 105°            |            |            | 35-41      | 42-46      | 47-50      | 51-52      |
| Lift activ                                        |            |            | 30-38      | 39-43      | 44-47      | 48         |
| Lift solid mechanisch und elektrisch              | 35         |            | 36-43      | 44-48      | 49-50      |            |
| Speedy 4teen L                                    | 26         |            | 27-35      | 36-40      | 41-44      | 45-51      |
| Speedy 4teen XL                                   | 26-27      |            | 28-35      | 36-40      | 41-45      | 46-52      |
| Speedy 4teen XXL                                  | 30         |            | 31-39      | 40-44      | 45-48      | 49-52      |
| Buddy classic D4-3                                | 27-33      | 34-42      |            |            |            |            |
| Buddy 4all D7                                     | 27-33      | 34-42      |            |            |            |            |
|                                                   |            |            |            |            |            |            |
| <b>Rohrlänge in mm</b>                            | <b>210</b> | <b>270</b> | <b>280</b> | <b>330</b> | <b>380</b> | <b>420</b> |

